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

Art. no. 983 003

GB

1. Reduces aching in your back, neck, shoulders, hips and joints
2. Increases blood circulation
3. Increases energy level
4. Gives a deep sensation of relaxation
5. Reduces muscle tension
6. Releases endorphins
7. Increases the bloods ability to transport oxygen
8. Improves sleep quality
9. Improves digestion
10. Improves mental health

How long?

Abilica Spike mat can be used as long as one wishes. There is no risk involved in e.g. sleeping on the mat, and children may use the mat if they wish to do so. For effective improvement of the blood circulation it is recommended to use the mat for at least 10 minutes. For people who have used the mat over several weeks, the normal time would be somewhere between 20 and 40 minutes for each use.

How often?

In general the mat may be used whenever one feels the need. The first period of use however, it is recommended to give the body a few days rest after using the mat for the first couple of weeks. This is to prevent the endorphins released when using the mat to gradually lose their effect.

How?

For people with extra sensitive skin it is recommended to use a T-shirt or a light sweater the first times they use the mat. The best result is still achieved on bare skin. Gently lie down on the mat and relax. Try to keep still, as excessive movement may cause the spikes to scratch the skin. The spikes should just act as pressure points.

Where?

Put a rolled up towel or small pillow under one end of the mat to support the head. It is also possible to place cushions under other parts of the mat to distribute the force from the spikes more evenly. Or simply place the mat on a softer surface, such as e.g. a bed. Placing the mat on a soft surface will result in a more gentle session.

- **IMPORTANT INFORMATION**

People suffering from haemophilia, or using blood diluting medicines, should not use the mat without first consulting a doctor. The pikes will not penetrate the skin, but minor internal muscular bleedings may occur.

- Do not use the mat during pregnancy.
- People suffering from psoriasis or severe eczema should consult a dermatologist before using the mat.
- The mat may be hand washed. Do not tumble dry.
- Do not fold or roll for storage.

SPORT SUPPLY®

Supplier:
Sport Supply Int. AS
P.O. Box 244
N-3051 Mjoendalen
Norway
E-post: post@sportsupply.no
www.abilica.com